# SPOT Software 中文操作


岑祥股份有限公司  
TSENG HSIANG LIFE SCIENCE LTD.

LIFE SCIENCE  
TECHNOLOGIES

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The image shows the splash screen for SPOT Software Version 4.6. The background is blue with a grid pattern and binary code. The text 'DIAGNOSTIC' is visible in the top left. The main text reads 'VERSION 4.6' in red and '4.6' in large green numbers. A yellow circular logo with 'SPOT' and a red crosshair is on the left. Below the logo, it says 'SPOTS SOFTWARE:'. Copyright information is listed at the bottom, including 'Copyright © 1997-2006 • Diagnostic Instruments, Inc.' and 'All rights reserved'. Technical details are provided in a table at the bottom.

Software:	Camera:
Version: 4.6.4.6/4.7.5.2	SPOT Insight 2Mp Mosaic 18.2 (247807)
Build Date: Feb 21 2007/Jan 17 2008	Serial Number: 247807
	Model: 18.2
	Revision: 0.0.1/0.0.19E

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The screenshot shows the SPOT software interface. The title bar reads "SPOT" and the menu bar includes "File", "Database", "View", "Camera", "Setup", "Macro", "Peripheral Devices", and "Help". The main window displays a splash screen with the following text:

*直接點選軟體 Spot software*  
*Creative control without the complex interface*

**windows:VERSION 4.6**  
ONE SOURCE FOR YOUR IMAGING SOLUTION

**SPOT** SPOT SOFTWARE: 4.6

Copyright© 1997-2005 • Diagnostic Instruments, Inc.  
All rights reserved.  
Portions Copyright© 1988-1998 • Sam Leffler  
Copyright© 1991-1996 • Silicon Graphics, Inc.  
Copyright© 2001 • David Taubman, University of New South Wales  
This software is based in part on the work of the independent JPEG Group.

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The bottom of the window shows a taskbar with the "start" button, several application icons, and the "SPOT Flex Mono 15.0 (101)" taskbar button. The system tray on the right shows "Brightfield" and the time "10:17 AM".

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The screenshot shows the SPOT software interface. The title bar reads "SPOT" and the menu bar includes "File", "Database", "View", "Camera", "Setup", "Macro", "Peripheral Devices", and "Help". The main workspace is a large grey area. On the right side, there is a vertical toolbar with icons for "Live", a camera, a magnifying glass, and a refresh symbol. At the bottom, the Windows taskbar is visible with the "start" button, several application icons, and the system tray showing "SPOT Flex Mono 15.0 (101)", "Brightfield", and the time "3:27 PM".

Four yellow callout boxes with red borders and red arrows pointing to specific areas of the interface provide the following information:

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
The screenshot displays the SPOT software interface. The main window, titled "Live Image - Running...", shows a histological image of tissue. On the right side of this window, there is a vertical "Focus Assist Indicator" bar, which is a gradient from blue at the top to red at the bottom. A red arrow points from a text box to this bar. Below the image, there is a "Focus Gauge" checkbox which is checked. Another red arrow points from a text box to this checkbox. The software has a menu bar (File, Database, View, Camera, Setup, Macro, Peripheral Devices, Help) and a taskbar at the bottom showing the Windows start button, several application icons, and the system tray with the time 9:48 AM.

SPOT

File Database View Camera Setup Macro Peripheral Devices Help

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
Brightfield

start Microsoft Power... SPOT 9:48 AM

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*Image setups let you quickly change applications*

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16 Megapixel  
2xbinBGSUB  
36 Megapixel  
4Megapixel  
64 Megapixel  
Brightfield  
CFP  
DAPI  
Darkfield  
DIC  
Factory Defaults  
FITC  
Fluorescence  
GFP  
Mono  
Rhodamine  
YFP  
Brightfield

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The screenshot displays the SPOT software interface. The main window shows a microscope image of a red, oval-shaped specimen on a black background. The software has a menu bar with options: File, Database, Edit, View, Camera, Setup, Macro, Peripheral Devices, Window, and Help. A toolbar on the right side includes icons for Live, image capture, zoom, and other functions. Two yellow callout boxes with red borders and arrows point to specific icons in the toolbar:

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The Windows taskbar at the bottom shows the Start button, several open applications (Microsoft PowerPoint, SPOT), and the system tray with the time 10:01 AM.

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The image shows the splash screen for SPOT Software Version 4.6. The background is blue with a grid pattern and binary code. The text 'DIAGNOSTIC' is in the top left, and 'VERSION 4.6' is in red at the top right. Below that, it says 'ONE SOURCE FOR YOUR IMAGING SOLUTION'. The SPOT logo is a yellow circle with 'SPOT' in black and red. To the right, '4.6' is written in large green numbers. Below the logo and version number, there is a copyright notice: 'Copyright © 1997-2006 • Diagnostic Instruments, Inc. All rights reserved. Portions Copyright © 1988-1996 • Sam Leffer Copyright © 1991-1996 • Silicon Graphics, Inc. Copyright © 2001 • David Taubman, University of New South Wales. This software is based in part on the work of the independent JPEG Group'. At the bottom, there is a table with software and camera information.

Software:	Camera: SPOT Insight 2Mp Mosaic 18.2 (247807)
Version: 4.6.4.6/4.7.5.2	Serial Number: 247807
Build Date: Feb 21 2007/Jan 17 2008	Model: 18.2
	Revision: 0.0.1/0.0.19E

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
The image shows a software splash screen for SPOT Software Version 4.6. The background is blue with a grid pattern and binary code. The text 'DIAGNOSTIC' is visible in the top left. The main title 'VERSION 4.6' is in red at the top right. Below it, 'ONE SOURCE FOR YOUR IMAGING SOLUTION' is written in white. A large yellow circular logo with 'SPOT' in black and red is on the left. To its right, 'SPOT SOFTWARE:' is written in white, followed by a large green '4.6'. A small yellow circular seal is also present. At the bottom, there is a white box containing software and camera information.

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The screenshot displays the SPOT software interface. The main window is titled "Live Image - Running...". On the left side of this window, there is a vertical menu with buttons for "Restart", "Pause", "Transfer", "Controls...", and "Histogram...". The "Controls..." button is circled in yellow. A yellow arrow points from this button to the "Live Image Controls" panel on the right. The "Live Image Controls" panel contains several sections: "Exposure" (Type: Auto/Manual, Binning: <none>, Exposure Time: 445.7916 msec), "Area" (Imaging Area: Full Chip/Selection, Region (from image setup), Always Display Full Chip Area), and "Post-Processing" (Gamma: 1.00, Display Images in Grayscale, Indicate Saturated Pixels). A red arrow points from the bottom of the "Live Image Controls" panel to a yellow text box. The text box contains the following text:

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The bottom of the screenshot shows the Windows taskbar with the "start" button, several application icons, and the system tray showing "SPOT Flex Mono 15.0 (101)", "Brightfield", and the time "8:55 AM".

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
The screenshot displays the SPOT software interface. The main window, titled "Live Image - Running...", shows a red, oval-shaped biological specimen against a black background. A yellow circle highlights the "Histogram..." button in the left-hand control panel. A yellow arrow points from this button to a smaller window titled "Live Image Histogram". This histogram window shows a graph with a peak on the right side, representing the intensity distribution of the image. A red arrow points from the histogram window to a yellow text box. The text box contains the following text:

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
The SPOT software interface includes a menu bar with options: File, Database, View, Camera, Setup, Macro, Peripheral Devices, and Help. The taskbar at the bottom shows the Windows Start button, several application icons, and the system tray with the time 9:58 AM.

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

The splash screen features a blue background with binary code and a large yellow '4.6' in the center. The SPOT logo is on the left, and the text 'SPOTS SOFTWARE: 4.6' is on the right. A gold seal is in the top right corner.

DIAGNOSTIC  
VERSION 4.6  
ONE SOURCE FOR YOUR IMAGING SOLUTION

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File Database Edit View Camera Setup Macro Peripheral Devices Window Help

Undo  
Copy to Clipboard Ctrl+C  
Resize...  
Rotate Right  
Rotate Left  
Flip Horizontal  
Flip Vertical  
Crop...  
Adjust RGB  
Adjust HSL  
Adjust HSV  
Filter  
Make Negative  
Convert to Gray Scale  
Match Color...  
Change Bit Depth...  
Subtract Image...  
Merge Images...  
Annotate  
Set Calibration...  
Add/Edit Calibration Marks...  
Add/Edit Measurements...  
Merge  
Memo...

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**Comprehensive Image Editing Tools**  
影像修飾工具列

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
The screenshot displays the SPOT software interface. The main window shows a fluorescence microscopy image of a cell with blue filaments and a red nucleus. A white line indicates a measurement of 67.55  $\mu\text{m}$ . A scale bar at the bottom right of the image indicates 20  $\mu\text{m}$ . The image is annotated with the following text:

Length: 67.55  $\mu\text{m}$

Texas Red - Histones (Ms)  
Alea Fluor 488 - PMP 70 (Rb)  
Alexa Fluor 360 - Phalloidin

20  $\mu\text{m}$

The software interface includes a menu bar (File, Database, Edit, View, Camera, Setup, Macro, Peripheral Devices, Window, Help) and a toolbar on the right with icons for Live, zoom, pan, and other functions. The Windows taskbar at the bottom shows the start button, SPOT application, and system tray with the time 9:18 AM.

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
Measurements 測量功能

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The screenshot displays the SPOT software interface. The main window shows a fluorescence microscopy image of a cell with a blue nucleus, red cytoplasm, and green filaments. A white line indicates a length of 67.55  $\mu\text{m}$ . A 20  $\mu\text{m}$  scale bar is present. The image is annotated with a yellow SPOT logo and a list of channels: Texas Red - Histones (Ms), Alexa Fluor 488 - PMP 70 (Rb), and Alexa Fluor 360 - Phalloidin. An inset image in the top right shows a zoomed-in view of the nucleus. The software's menu bar includes File, Database, Edit, View, Camera, Setup, Macro, Peripheral Devices, Window, and Help. The status bar at the bottom shows 'SPOT Flex Mono 15.0 (101)', '1005x802 24bpp', '(982,470) 7.0,105', 'Fit to W', 'Modifi', and 'Rhodamine'. The Windows taskbar at the bottom shows the start button, several icons, and the time 9:18 AM.

Length: 67.55  $\mu\text{m}$

20  $\mu\text{m}$

Texas Red - Histones (Ms)  
Alexa Fluor 488 - PMP 70 (Rb)  
Alexa Fluor 360 - Phalloidin

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
The screenshot displays the SPOT software interface with three image windows open. The top-left window shows a red oval-shaped object on a black background. The top-right window shows a green speckled pattern on a black background. The bottom-left window shows a blue striped pattern on a black background. The software interface includes a menu bar (File, Database, Edit, View, Camera, Setup, Macro, Peripheral Devices, Window, Help) and a toolbar on the right side. The taskbar at the bottom shows the Windows Start button, several application icons, and the system tray with the time 9:10 AM.

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
The image shows a software splash screen for SPOT Software Version 4.6. The background is blue with a grid pattern and binary code. The text 'DIAGNOSTIC' is in the top left, and 'VERSION 4.6' is in red at the top right. Below that, it says 'ONE SOURCE FOR YOUR IMAGING SOLUTION'. The SPOT logo is a yellow circle with a red crosshair. To its right, 'SPOT SOFTWARE:' is written in white, followed by a large '4.6' in green. A small circular seal is also visible. At the bottom, there is a white box with software and camera information.

DIAGNOSTIC

VERSION 4.6

ONE SOURCE FOR YOUR IMAGING SOLUTION

SPOT


SPOT SOFTWARE:

4.6

Copyright © 1997-2006 • Diagnostic Instruments, Inc.  
All rights reserved.  
Portions Copyright © 1988-1996 • Sam Leffer  
Copyright © 1991-1996 • Silicon Graphics, Inc.  
Copyright © 2001 • David Taubman, University of New South Wales  
This software is based in part on the work of the independent JPEG Group

Software:	Camera: SPOT Insight 2Mp Mosaic 18.2 (247807)
Version: 4.6.4.6/4.7.5.2	Serial Number: 247807
Build Date: Feb 21 2007/Jan 17 2008	Model: 18.2
	Revision: 0.0.1/0.0.19E

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The screenshot displays the SPOT software interface. The main window shows a large fluorescence image of a rat brain section with the following text overlay: "Alexa Fluor 568 - GFAP (Ms)", "Alexa Fluor 488 - MAP-2 (Ch)", "Hoechst 33342", and "Rat Brain Horizontal 8-Micrometer Section". Below this image, there is a yellow text box with the text: "Search, Open, Copy, Delete and Print images quickly and easily" and "可任意搜尋, 開啓, 複製, 刪除, 列印等相關功能". To the right of the main image is a "File Images" window showing a grid of 46 thumbnail images, labeled from "ratbrain\_largeAL-SMI-312.jpg" to "ratcoloa\_large.jpg". The software interface includes a menu bar (File, Database, View, Camera, Setup, Macro, Peripheral Devices, Window, Help) and a taskbar at the bottom with the Windows Start button, taskbar icons for "Microsoft Power...", "Imaging", and "SPOT", and a system tray showing "Brightfield" and the time "2:55 PM".

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
The screenshot displays the SPOT software interface. On the left, a large window shows a fluorescence microscopy image of a cell with blue filaments and a red nucleus. On the right, a 'Print - HP Photosmart 8400 Series' dialog box is open, showing a preview of three images stacked vertically. The dialog includes a 'Page: 1 of 1' indicator, a 'Print' button, and a 'Setup...' button. Below the preview, the 'Printable Area' is listed as 7.992 x 10.315. The 'Top Image' section has input fields for 'Left: 2.194', 'Width: 4.077', 'Top: 6.995', and 'Height: 3.438'. There are also checkboxes for 'Annotations', 'Calibration Marks', and 'Measurements'. The 'Show:' section has 'Annotations' and 'Calibration Marks' checked. The 'Add...', 'Remove', and 'Arrange...' buttons are at the bottom of the dialog. The status bar at the bottom shows 'SPOT Flex Mono 15.0 (101)', '1005x802 24bpp', '(1003.639) 0,0,2', 'Fit to W', and 'Rhodamine'.

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The screenshot displays the SPOT software splash screen. The window title is "SPOT" and the menu bar includes "File", "Database", "View", "Camera", "Setup", "Macro", "Peripheral Devices", and "Help". The main content area features the following text:

**Creative control without the complex interface**

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**windows:VERSION 4.6**

ONE SOURCE FOR YOUR IMAGING SOLUTION

**SPOT** SPOT SOFTWARE: 4.6

Copyright © 1997-2005 • Diagnostic Instruments, Inc.  
All rights reserved.  
Portions Copyright: 1995-1996 • Sam Laffin;  
Copyright: 1993-1994 • Silicon Graphics, Inc.  
Copyright: 2001 • David Tuckman, University of New South Wales.  
This software is based in part on the work of the independent JPECO Group.

Software:	Camera:
Version: 4.6.0.1 Prelim/4.6.0.1	Serial Number: 101
Build Date: Oct 7 2005/Oct 7 2005	Model: SPOT Flex Mono
	Revision: 0.0.1/0.0.26 Eng


At the bottom of the splash screen, it says "Spot software included with all Diagnostic Instruments Cameras". The Windows taskbar at the bottom shows the "start" button, several icons, and the system tray with the time "10:17 AM".

*Spot software included with all Diagnostic Instruments Cameras*

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Get Sequential Images

Number of Images: 3

Compute Exposure before Each Image

Time Interval:  As Short as Possible  
0 hour 0 min 2.000 sec

Place Images in an Image Sequence

Defer Image Processing/Display until Done


Status  
Images Acquired: 0      Processed: 0  
Waiting to begin...

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
The screenshot displays the SPOT software interface. A window titled "new image sequence 1 (3/3)" shows a video frame of a document. Overlaid on this is the "Get Sequential Images" dialog box. The dialog box contains the following settings:

- Number of Images: 3
- Compute Exposure before Each Image
- Time Interval:  As Short as Possible (0 hour, 0 min, 0.000 sec)
- Place Images in an Image Sequence
- Defer Image Processing/Display until Done
- Status: Images Acquired: 0, Processed: 0, Waiting to begin...

Buttons for "Begin", "Close", and "Help" are visible at the bottom of the dialog box. The software's status bar at the bottom indicates "SPOT Insight 2Mp Mosaic 18.2 (247807)", resolution "1600x1200 24bpp", coordinates "(712,49) 11,134,92", and a "Preview" dropdown menu. The system tray shows the time as "下午 03:18".

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
The screenshot displays the SPOT software interface. The main window shows a video sequence titled "new image sequence" with a context menu open over it. The menu options are:

- Extract Images...
- Insert Image...
- Append Image...
- Remove Images...
- Move Images...
- Play Forward (Ctrl+Space)
- Play Reverse (Ctrl+Backspace)
- Stop Play (Enter)
- Step Forward (Space)
- Step Reverse (Backspace)

The video content shows a scene with yellow and white papers, possibly a document or a book. The interface includes a menu bar (File, Database, Edit, Sequence, View, Camera, Setup, Macro, Peripheral Devices, Window, Help) and a toolbar on the right side. The taskbar at the bottom shows the SPOT application and Microsoft PowerPoint.

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此影像序列檔New Images sequence, 也可以放映操作, 點選View/  
Show Sequence Control功能鍵

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The screenshot displays the SPOT software interface. The main window shows a video preview of a cluttered desk with papers and a yellow envelope. An "Image Sequence Controls" dialog box is open, featuring playback navigation buttons (stop, previous, play, next, first, last), an "Index Number" section with "First: 1", "Current: 3", and "Last: 3", and "Play Options" including "Interval: 150" (msec/sec), "At End: Stop/Restart/Reverse", and "Close/Apply/Help" buttons. A vertical toolbar on the right includes a "Live" button and various tool icons. The Windows taskbar at the bottom shows the SPOT application, a Microsoft PowerPoint window, and the system clock at 03:19.

new image sequence 1 (3/3)

Image Sequence Controls

Index Number  
First: 1 Current: 3 Last: 3

Play Options  
Interval: 150 msec sec  
At End: Stop Restart Reverse  
Close Apply Help

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
new i...

SPOT Insight 2Mp Mosaic 18.2 (247807) 1600x1200 24bpp (1490,365) 2,134,96 Fit to Win View 預覽

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The screenshot shows the SPOT software interface. The 'File' menu is open, and 'Export Movie File...' is highlighted. The main window displays a video frame of a desk with papers and a yellow envelope. The sidebar on the right contains various tool icons. A text box on the right explains the 'Export Movie File' option.

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