# **LightCycler<sup>®</sup> 480 Software v1.5.0**

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Pos	Color	Repl Of	Sample Name
A1	Blue		Sample 1
A2	Red		Sample 1
A3	Green		Sample 3
A4	Magenta		Sample 4
A5	Grey		Sample 5
A6	Yellow		Sample 6
A7	Brown		Sample 7
A8	Cyan	A8	Treat
A9	Dark Green	A8	Treat
A10	Orange		Sample 10
A11	Purple		Sample 11
A12	Olive		Sample 12

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Pos	Color	Repl Of	Sample Name
A1	Blue		Sample 1
A2	Red		Sample 1
A3	Green		Sample 3

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絕對定量 (Abs Quant) :

Pos	Color	Repl Of	Sample Name	Quantification Sample Type	Concentration
A1	Blue		Sample 1	Unknown	
A2	Red		Sample 1	Unknown	
A3	Green		Sample 3	Unknown	
A4	Magenta		Sample 4	Unknown	
A5	Grey		Sample 5	Positive Control/Calibrator	
A6	Yellow		Sample 6	Negative Control Standard	

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相對定量 (Rel Quant) :

Pos	Color	Repl Of	Sample Name	Combined Sample and Target Type	Concentration	Target Name	Efficiency
A1	Blue		Sample 1	Unassigned Unknown			2.00
A2	Red		Sample 2	Unassigned Unknown			2.00
A3	Green		Sample 3	Unassigned Unknown			2.00
A4	Magenta		Sample 4	Unassigned Unknown			2.00
A5	Grey		Sample 5	Target Unknown			2.00
A6	Yellow		Sample 6	Target PosCalibrator			2.00
A7	Brown		Sample 7	Target Negative			2.00
B1	Blue		Sample 13	Target Standard			2.00
B2	Red		Sample 14	Ref Unknown			2.00
B3	Green		Sample 15	Ref PosCalibrator			2.00
B4	Magenta		Sample 16	Ref Negative			2.00
B5	Grey		Sample 17	Ref Standard			2.00
B6	Yellow		Sample 18	Unassigned Unknown			2.00
B7	Brown		Sample 19	Unassigned PosCal			2.00
				Unassigned Neg			2.00
				Unassigned Std			2.00

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
The screenshot shows the software interface for editing a sample. It is divided into three main sections:


- Step 1: Select Workflow:** Contains radio buttons for 'Abs Quant', 'Rel Quant' (selected), 'Scanning', 'Color Comp', 'Tm', 'Melt Geno', and 'Endpt Geno'. A label '分析模式' (Analysis Mode) points to this section.
- Step 2: Select Samples:** Features a 'Subset' dropdown set to 'All Samples', a grid of 12 columns and 6 rows (A-F) with 'U' in each cell, and various icons for selection. A label '選取樣品' (Select Sample) points to this section.
- Step 3: Edit Rel Quant Properties:** Includes fields for 'Sample Name', 'Sample Type' (Unknown, Negative Control, Positive Control/Calibrator, Standard), 'Gene target' (Target, Reference, Unassigned), and 'Eff' (2.00). A 'Make Replicates' button is at the bottom. A label '編輯樣品名稱與定義樣品' (Edit sample name and define sample) points to the Sample Name and Type fields, and another label '編輯基因名稱與 PCR 效率' (Edit gene name and PCR efficiency) points to the Gene target and Eff fields.

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**Abs Quant / 2nd Derivative Max**

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**Abs Quant / Fit Points**

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**Advanced Relative Quantification**

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Results			Target Name
Target Name	Filter Combination	Standards Efficiency	Efficiency Value
Target 1	465-510	Efficiency	1.80
Target 2	465-510	Efficiency	2.00
Reference 1	465-510	Efficiency	2.00
Reference 2	465-510	Efficiency	2.00

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**Ratios**

- Display Target/Reference Ratio
- Display Normalized Ratio
- Show Ratio Errors

If checked, the ratio shows in the result table and the bar chart.

**Result Table**

- Display Median/Mean Cp
- Display Cp Error

Checked values are displayed as columns in the result table.

**Manual Pair Editing**

- Enable Manual Pair Editing

If checked, the Manual Pair Editing Tab is visible.

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

The screenshot shows the software interface with the following elements:

- Report Settings:** A list of items to include in the report, with a red box around the list and the 'Default Settings' button below it.
 

Item	Checked
Experiment	<input checked="" type="checkbox"/>
Protocol	<input checked="" type="checkbox"/>
Samples	<input checked="" type="checkbox"/>
Instrument	<input checked="" type="checkbox"/>
Revision History	<input checked="" type="checkbox"/>
Abs Quant/2nd Derivative Me	<input checked="" type="checkbox"/>
Settings	<input checked="" type="checkbox"/>
Results	<input type="checkbox"/>
Statistics	<input type="checkbox"/>
Amplification Curves	<input type="checkbox"/>
Standard Curve	<input type="checkbox"/>
- View Report:** A window showing report details and a table of parameters. The 'PDF' button is circled in red.
 

Experiment	
Creation Date	2008/1/30 下午 03:11:11
Last Modified Date	2008/2/5 上午 10:30:37
Operator	System Admin
Owner	System Admin
Start Time	2008/1/30 下午 03:51:03
End Time	2008/1/30 下午 04:35:51
Run State	Completed
Software Version	LCS480 1.2.0.169
Macro	Macro Owner
Macro Status	
Templates	Roche reagent384
Plate ID	
Run Notes	

Programs							
Program Name pre-incubation							
Cycles	1	Analysis Mode	None				
Target (°C)	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)	Acquisition Mode
95	00:05:00	4.8		0	0	0	None

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