

【尿液收集方法】

[Guidelines for collection of Urine Specimens]

【尿液常規檢驗】:

- 1. 請先排掉前段尿液。
- 2. 用紙尿杯接取中段尿液。
- 3. 將紙杯尿液倒入紅蓋尿液收集管中。 尿液常規檢查約八分滿(12毫升)。

【生化尿液檢驗】:

- 1. 請先排掉前段尿液。
- 2. 用紙尿杯接取中段尿液。
- 3. 將紙杯尿液倒入白蓋尿液收集管中。生化尿液檢驗約試管一半的量。

【尿液細菌培養】:

- 1. 採集尿液前,請先以優碘棉片清潔尿道口。
- 2. 請先排掉前段尿液。
- 3. 以無菌尿杯直接接取 20-30 毫升中段尿液。

不可先排尿於紙尿杯後再倒入無菌尿杯。

[Urine Routine Collection]:

- 1. Please drain the front part of the urine first.
- 2. Use the paper urine cup to catch the middle part of the urine.
- 3. Pour the paper cup urine into the red cap urine collection tube. The urine should be about 80% full (10-12 ml) for routine examination.

[Urine Biochemistry Test Collection]:

- 1. Please drain the front part of the urine first.
- 2. Use the paper urine cup to collect the middle part of the urine.
- 3. Pour the paper cup urine into the white cap urine collection tube. The urine should be about <u>half</u> full for Urine Biochemistry Tests.


(Urine Culture Collection):


- 1. Before collecting urine, clean the urethral meatus with Povidone-Iodine Prep Pad.
- 2. Please drain the front part of urine first.
- 3. Use the sterile urine cup to receive 20-30 mL of middle part urine directly.


DO NOT urinate into the paper cup first then pour urine into the sterile container.

◎注意事項:

- 1. 女性正值月經期間,原則上應避免做尿液檢驗。
- 2. 尿液檢體請在1小時內送達檢驗室。
- 3. 如無法立刻送檢,應將檢體置於冷藏環境(2-8℃),且最好不要超過4小時
- 4. 請您攜帶健保卡、申請單與檢體,於服務時間送至收檢櫃檯。
- 5. 若有任何疑問,請來電詢問臨床病理科。
- 1. Female should avoid urine tests during menstruating period.
- 2. Urine test should be sent to the laboratory within 1 hour.
- 3. If you are unable to send the test immediately, the sample should be kept in the refrigerator for less than 4 hours.
- 4. Please bring your National Health Insurance Card, application form and sample to the laboratory during service hours.
- 5. If you have further queries, please do not hesitate to contact us. Thank you.







Neihu District: (02)87923311

#17438 or 88094

汀州院區:(02) 23659055 轉分機 88617 Ting-jhou District: (02)23659055 #88617

